



Yes, you can.™

## **➤ INVACARE® PRONTO M5 I**

Presented in a stylish compact size, the Invacare Pronto M5 I power wheelchair with Surestep has true centre-wheel drive performance for exceptional manoeuvrability and intuitive driving. Featuring a light grey semi-recline van seat and innovative Surestep technology, the Pronto M5 I provides a smooth, less jarring ride over transitions and thresholds up to 5cm in height while maintaining stability. The Pronto M5 I is attractively designed for optimum manoeuvrability and accessibility.

*The M5 I is an affordable, centre wheel drive packed with the latest technology.*



## ► SureStep® Technology

SureStep technology provides optimum traction and stability as the Pronto M51 drives over transitions and uneven ground up to two inches in height. The powerful motors lift the front castors to assist with changes in height for a smooth, less jarring ride over rough surfaces. SureStep technology allows six wheels to remain in contact with the ground to maximize performance, improve stability and enhance comfort for a more controllable power chair.

## ► Centre Wheel Drive


True centre wheel drive provides the most intuitive driving, an astounding turning radius and exceptional manoeuvrability.

## ► Stylish New Semi-Recline Van Seat


Generous new contours provide support where needed, seat reclines for comfort, and backrest extends to fit any size consumer. Seat to floor as low as 48 cm with new Light Grey Semi-Recline van seat - provides ease of transferring into and out of the chair! Seat to floor is adjustable up to 58cm.



### TECHNICAL DATA:




61 cm/24" (inches)  
(excludes joystick)



41-51 cm/  
16-20" (inches)




41-46 cm/  
16-18" (inches)



107 cm/42" (inches)  
approx. 115 cm/45" (inches)  
with the headrest extended



86 cm/34" (inches)  
with footboard up  
97 cm/38" (inches)  
with footboard  
down




90 kg / 198 lbs with  
batteries  
69 kg / 151 lbs without  
batteries



4 mph (6.4 kph)  
(max)



136 kg / 300 lbs



50 cm / 19.5" (inches)  
front & rear turning



9° (degrees)

### Limited Warranties

Base frame	60 months
Seat frame	60 months
Electronics	12 months
Motor/gearbox	12 months

### FEATURES:



**Front castors**  
lift up to accommodate  
transitions up to two-inches  
for a smooth and  
controllable ride.



**Optional Extra:  
Rear Bag Kit**  
For extra storage.

### FRAME COLOUR:




Ruby Metallic Red



Deep Blue

### UPHOLSTERY:



Light Grey

Pronto M51-NZ/AU-05/06

Invacare® International Sarr 2005.  
© All rights reserved.

# INVACARE

Yes, you can.™

### Local Distributor:

Invacare Australia Pty Ltd ABN 45 074 676 378.  
1 Lenton Place, PO Box 5002, North Rocks, NSW 2151, Australia.  
Freephone: 1 800 460 460, Fax: 02 8839 5353,  
Email: sales@invacare.com.au, www.invacare.com.au

Invacare New Zealand,  
4 Westfield Place, PO Box 62-124, Mt Wellington, Auckland, New Zealand.  
Freephone: 0508 468 222, Freefax: 0508 807 788  
Email: sales@invacare.co.nz, www.invacare.co.nz

Invacare Asia Ltd,  
1 Lenton Place, PO Box 5002, North Rocks, NSW 2151, Australia.  
Phone: +61 2 8839 5330, Fax +61 2 8839 5331,  
Email: asiasales@invacare.com, www.invacare.com

All information is subject to change without prior notice. Photographs are for illustration purposes only and may vary from product description. Invacare® is a registered trademark and Yes, you can.™ is a trademark of Invacare® Corporation. 2005